


The arching dimension values in the attached report indicate that the hopper opening must be greater than 170mm for Milk Powder 6 in order to ensure reliable discharge of all powders. This dimension also applies to any conveyance device directly downstream of the hopper, such as the spacing between flutes on a screw feeder.

Data from the shear cell can also be used to calculate the critical rathole dimension, but this requires user input on bin shape, cross section dimension(s), and hopper type (wedge, conical). With the assumption of a 2m bin diameter, the rathole diameter was calculated for Milk Powder 6 at 575mm. The other milk powders all had lower values for rathole diameter.

The milk powder industry has become one of the major users of shear cells. Undoubtedly the test information has proven more useful in solving the problems of predicting flowability and designing the processing equipment, such as hopper opening and half angle, to handle the variations in raw materials and final blend requirements.

Bulk Density Graph (230cc Trough)



This bulk density graph illustrates the density of the samples over different consolidating stresses. All 6 milk powders show a moderate to high degree of compressibility. Loose fill density values start between 0.38 and 0.41 g/cm³. Milk powder 3 shows the least change, approximately 45%, while Milk Powder 6 shows the greatest change, approximately 70%. Powders with the greatest amount of change in density are more likely to exhibit flow behavior problems.

Author: Robert G. McGregor, General Manager -Global Marketing

Brookfield Engineering Laboratories, Inc., Middleboro, MA 02346

1.946.6200 Ext. 7143 Email: r_mcgregor@brookfieldengineering.com

Website: <http://www.brookfieldengineering.com>

Reprint Permission VertMarkets - Pharmaceutical Online - January 2011



VISCOMETERS

RHEOMETERS

TEXTURE ANALYZERS

POWDER FLOW TESTERS



T: 800.628.8319 / 508.946.6200
www.brookfieldengineering.com